

# Phonak Audéo™ 1.

## Product information

Phonak Audéo I is the comprehensive RIC portfolio for mild to profound hearing losses. Powered by the innovative technology of the Phonak Infinio platform, Audéo I hearing aids are the solution for clients that prefer flexible and discreet RIC housings.



# Product description



**Audéo I-Sphere**  
Audéo I-Sphere Trial  
190/170



**Audéo I-R**  
Audéo I-R Trial  
190/170/150/130

Battery		built-in rechargeable	built-in rechargeable
Multi-function button		•	•
Bluetooth®		•	•
Hydrophobic coating		•	•
IP rating		IP68 <sup>1</sup>	IP68 <sup>1</sup>
Dimensions	L x W x D	31.9 x 14.0 x 9.5 mm (1.26 x 0.55 x 0.37")	28.8 x 12.4 x 8.7 mm (1.13 x 0.49 x 0.34")
	Weight	3.0 g (0.11 oz)	2.4 g (0.08 oz)
Charging times <sup>2</sup>	0% to 100%	3.5 h	3 h
	0% to 70%	2 h	2 h
	0% to 40%	1 h	1 h
	30 min. fitting	15 min.	15 min.
<b>S Receiver 6.0</b>			
Max. Power Output (dB SPL)	2 cc coupler	112	112
Max. gain (dB)	2 cc coupler	47	47
Frequency range (Hz)		<100 - 10000	<100 - 10000
<b>M Receiver 6.0</b>			
Maximum Power Output (dB SPL)	2 cc coupler	114	114
Max. gain (dB)	2 cc coupler	51	51
Frequency range (Hz)		<100 - 10000	<100 - 10000
<b>MAV Receiver 6.0</b>			
Maximum Power Output (dB SPL)	2 cc coupler	119	119
Max. gain (dB)	2 cc coupler	50	50
Frequency range (Hz)		<100 - 9000	<100 - 9000
<b>P Receiver 6.0</b>			
Maximum Power Output (dB SPL)	2 cc coupler	122	122
Max. gain (dB)	2 cc coupler	59	59
Frequency range (Hz)		<100 - 7500	<100 - 7500
<b>UP Receiver 6.0</b>			
Max. Power Output (dB SPL)	2 cc coupler	132	132
Max. gain (dB)	2 cc coupler	71	71
Frequency range (Hz)		<100 - 8000	<100 - 8000

## Housing colors



**P1**  
Sand  
Beige



**P3**  
Sandalwood



**P4**  
Chestnut



**P5**  
Champagne



**P6**  
Silver  
Gray



**P7**  
Graphite  
Gray



**P8**  
Velvet  
Black



**1A**  
Copper



**8A**  
Sienna  
Red



**9A**  
Atlantic  
Green

<sup>1</sup> IP68 indicates that the hearing aid is water and dust resistant. It survived continuous immersion in 1 meter of fresh water for 60 minutes and 8 hours in a dust chamber as per the IEC60529 standard.

<sup>2</sup> Actual charging time depends on the remaining battery life.

The Bluetooth® word mark and logos are registered trademarks owned by Bluetooth SIG, Inc. and any use of such marks by Sonova AG is under license.

# Performance levels

AutoSense OS™ 6.0	Audéo I-Sphere		Audéo I-R			
	Premium	Advanced	Premium	Advanced	Standard	Essential
Spheric Speech in Loud Noise	•	•				
Speech in Loud Noise			•			
Music	•	•	•	•		
Speech in Car	•	•	•	•	•	
Comfort in Noise	•	•	•	•	•	
Speech in Noise	•	•	•	•	•	•
Calm Situation	•	•	•	•	•	•
Media speech + mic	•	•	•	•	•	•
Media music + mic	•	•	•	•	•	•
<b>Additional programs</b>						
Max. additional programs	4	4	4	4	4	4
Spheric Speech in Loud Noise	•	•				
Speech in Loud Noise	•	•	•	•		
Speech in 360°	•	•	•	•		
Speech in Car	•	•	•	•	•	
Speech in Noise	•	•	•	•	•	•
Calm Situation	•	•	•	•	•	•
Comfort in Noise	•	•	•	•	•	•
Music	•	•	•	•	•	•
Acoustic phone	•	•	•	•	•	•
Custom program	•	•	•	•	•	•
<b>Streaming programs</b>						
RogerDirect + mic	•	•	•	•	•	•
PartnerMic + mic	•	•	•	•	•	•
Phone call + mic	•	•	•	•	•	•
<b>Features</b>						
Spheric Speech Clarity	•	•				
Speech Enhancer	•		•			
SpeechSensor	•	•	•			
StereoZoom 2.0	•	•	•	•		
Dynamic Noise Cancellation	•	•	•	•		
Motion Sensor Hearing	•	•	•	•	•	
RogerDirect™	•	•	•	•	•	•
ActiveVent™	•	•	•	•	•	•
Tap control						
- Siri/Google Assistant	•	•	•	•		
- Pause/resume streaming	•	•	•	•		
- Accept/end phone call	•	•	•	•		
SoundRelax	•	•	•	•	•	
DuoPhone	•	•	•	•	•	
Real Ear Sound	•	•	•	•	•	•
UltraZoom	•	•	•	•	•	•
Health functionalities <sup>1</sup>	•	•	•	•	•	•
WindBlock	•	•	•	•	•	•
User Preference Tuning	•	•	•	•	•	•
Tinnitus Balance	•	•	•	•	•	•
Auto Acclimatization	•	•	•	•	•	•
Environmental Balance	•	•	•	•	•	•
SoundRecover2	•	•	•	•	•	•
WhistleBlock	•	•	•	•	•	•
NoiseBlock	•	•	•	•	•	•
QuickSync	•	•	•	•	•	•
AOV	•	•	•	•	•	•
Fine tuning channels	20	20	20	20	16	12

SmartSpeech Technology

<sup>1</sup> Available in combination with myPhonak

# CROS compatibility

	Audéo I-Sphere Audéo I-Sphere Trial	Audéo I-R Audéo I-R Trial
Phonak CROS I-R & Phonak CROS I-R Trial		•


## Fitting

	Audéo I-Sphere Audéo I-Sphere Trial	Audéo I-R Audéo I-R Trial
Phonak Target	10.0 or higher	10.0 or higher
Noahlink Wireless, Noahlink Wireless 2	•	•

## Acoustic coupling

Length	S Receiver	M Receiver	MAV Receiver	P Receiver	UP Receiver
Available in five different lengths (00, 0, 1, 2, 3), left and right	•	•	• <sup>1</sup>	•	• <sup>1</sup>
<b>Dome</b>					
Cap dome	•	•		•	
Open dome Available in three sizes (S, M, L)	•	•		•	
Vented dome Available in three sizes (S, M, L)	•	•		•	
Power dome Available in three sizes (S, M, L)	•	•		•	
<b>SlimTip</b>					
Acrylic	•	•		•	
Silicone	•	•		•	
Titanium AV			• <sup>1</sup>		
Universal AV Available in four shapes (A, B, C, D)			• <sup>1</sup>		
<b>cShell</b>					
Acrylic	• <sup>1</sup>	• <sup>1</sup>		• <sup>1</sup>	• <sup>1</sup>
Titanium	• <sup>1</sup>	• <sup>1</sup>		• <sup>1</sup>	• <sup>1</sup>

### Fitting ranges



Mild to moderately-severe hearing loss.

Mild to severe hearing loss.

Mild to severe hearing loss.

Mild to severe hearing loss.

Moderate to profound hearing loss.

<sup>1</sup> Length option 00 not available

# Phonak Audéo I charging options



Phonak Charger RIC I

Phonak ChargerGo RIC I

Phonak ChargerGo RIC Sphere I

Compatible with Audéo I-R and Audéo I-R Trial	•	•	
Compatible with Audéo I-Sphere and Audéo I-Sphere Trial			•
Charger with integrated battery		•	•

## Phonak Audéo I Bluetooth® version

Direct calls/audio streaming

Bluetooth 5.3 wireless technology and most older Bluetooth phones

## Wireless communication portfolio

	Audéo I-Sphere Audéo I-Sphere Trial 190/170	Audéo I-R Audéo I-R Trial 190/170/150/130
<b>Wireless accessory</b>		
TV Connector	•	•
Phonak PartnerMic	•	•
Phonak RemoteControl	•	•
<b>Roger</b>		
Roger microphone portfolio	•	•
<b>Apps</b>		
myPhonak	•	•



For more information, visit [www.phonak.com/audeo-infinio](http://www.phonak.com/audeo-infinio)

For information on Roger microphones, visit [www.phonak.com/roger](http://www.phonak.com/roger)

For more information on myPhonak and myPhonak Junior, visit [www.phonak.com/apps](http://www.phonak.com/apps)

# life is on

Since 1947, Phonak is dedicated to preserving life quality by opening new acoustic worlds. Already back then, in the days of our foundation, our company was driven by a formative conviction: We believe that well-hearing equates to well-being and thus is essential for living life to the fullest. In fact, the sense of hearing is directly linked to social, emotional, cognitive and physical well-being. Today as in future, we strive to offer the broadest portfolio of innovative hearing solutions. And, together with our hearing care professionals, we keep on focusing on what matters most: improving speech understanding, changing people's lives and having a positive effect on society as a whole.

[phonak.com](https://phonak.com)



027-0712-02/M1.00/2024-05/na © 2024 Sonova AG All rights reserved